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

HYALREPAIR®

HYALREPAIR®-06 Bioreparant

Gel based on modified hyaluronic acid hydrous sterile

Usage instruction

Ingredients

Sodium hyaluronate, copolymer of hyaluronic acid with sodium ascorbyl phosphate, copolymer of hyaluronic acid with riboflavin, sodium chloride, phosphate buffer, water for injections.

Description

HYALREPAIR®-06 Bioreparant is a transparent, yellow colour, homogeneous, gel-like substance supplied in a dark glass vial.

A unique formula of hyaluronic acid modification in form of copolymers with sodium ascorbyl phosphate and riboflavin. Contains high-purity hyaluronic acid obtained by bacterial fermentation.

Mechanism of action

Thanks to the fact that in HYALREPAIR®-06 Bioreparant hyaluronic acid modification is able to create copolymers with other substances, the period of its biodegradation slows down making hyaluronic acid not recognizable for hyaluronidases. This contributes to longterm retention in tissues of hyaluronic acid, vitamin C and riboflavin also, which are linked to it with solid phase modification. During biodegradation of HYALREPAIR®-06 Bioreparant in tissues hyaluronic acid, vitamin C and riboflavin are gradually released to restore the intercellular matrix and provide cellular activity (the production of proper young collagen, elastin and hyaluronic acid).

HYALREPAIR®-06 Bioreparant creates an optimum environment for normal functioning of fibroblasts which are responsible for synthesis and secretion of intercellular matrix proteins (collagen, elastin) and mucopolysaccharides, including their own hyaluronic acid.

Hyaluronic acid is a natural component of human skin, ensures the maintenance of three-dimensional organization of the components of the intercellular matrix and improves the process of natural metabolism, participates in transport and distribution of water in tissues, determines barrier and protective functions of intercellular space, improves skin elasticity, normalizes its moisture content. Hyaluronic acid

in form of the copolymers in HYALREPAIR®-06 Bioreparant effects targeted delivery of ascorbic acid and riboflavin to skin cells.

Gradual release of sodium ascorbyl phosphate from a copolymer of hyaluronic acid in HYALREPAIR®-06 Bioreparant maintains for a long time its physiological concentration in tissues. When the sodium ascorbyl phosphate is split, the released vitamin C (ascorbic acid) does participate in the synthesis of collagen and through lysine it participates in the formation of lysine bridges in the collagen structure. Similarly, ascorbic acid, converting lysine to oxylysine, participates in the formation of cross-linking of elastin and collagen and, thereby, stabilizes the reticular matrix of connective tissue. An important property of ascorbic acid is its unique ability to potentiate the formation of mucopolysaccharides of connective tissue (hyaluronic and chondroitin-sulfuric acids). This contributes to the formation of the support skin layer and the rapid tissue regeneration. Vitamin C protects skin keratinocytes from UV radiation, influences the process of melanin synthesis, thus decreasing hyperpigmentation. Ascorbic acid is a powerful antioxidant. Ascorbic acid regulates hydrogen ions transport in a number of biochemical reactions, improves utilization of glucose in tricarboxylic acid cycle, contributes to maintenance of colloidal state of intercellular substance and normal capillary permeability (due to inhibition of hyaluronidase activity). Vitamin C is necessary for protection against viral and bacterial infections, stimulation of macrophages, induction of endogenous interferon. Physiological doses of ascorbic acid stimulate phagocytosis.

Riboflavin stabilizes the intercellular matrix of connective tissue, it facilitates absorption of oxygen by skin cells and accelerates transformation of pyridoxine into its active form. Riboflavin participates in oxidation-reduction processes, activates synthesis of flavin coenzymes which accelerate tissue regeneration processes and interacts with adenosine triphosphate (ATP), the main carrier of energy in the cell.

Indications for use

1. Correction of age-related skin changes: restoration of lost volume of the dermis, filling of depressions and wrinkles.
2. Restoration of skin hydration.
3. Strengthening of blood vessels.
4. Hyperpigmentation.
5. Age-related decrease in skin elasticity in combination with tendency to onset of inflammatory elements.
6. Acne.
7. Change in micro and macro relief (enlarged pores, areas of hyperkeratosis).
8. Reactivity and sensitivity of the skin.

9. Atrophic scars (striae, post-acne).
10. Xerosis.
11. Seborrhea, including the scalp.
12. Diffuse hair loss.

Contraindications

General:

1. Acute inflammatory eruptions (herpes) or exacerbation of chronic skin diseases in supposed injection areas.
2. Autoimmune diseases.
3. Chronic diseases in exacerbation and decompensation stages.
4. Febrile state.
5. Hypersensitivity to the components.
6. Diseases involving impairment of healing processes.
7. Blood-clotting disorders (hemophilia).

Relative:

1. Epidermis injuries in the application area (e.g. post-peeling period).
2. Administration of anticoagulants or blood platelet aggregation inhibitors (e.g. ASS).
3. Morbid fear of needles.
4. Overstated patient expectations.

Precautions

Important! It is vital to comply with the injection guide to ensure due result of the procedure. Injection of HYALREPAIR®-06 Bioreparant may only be performed by specialists who have been duly authorized there to according to the local legislation and properly trained. Patients should be informed that the areas of injection should not be exposed to intensive heating (sun exposure, artificial tanning, sauna etc.) and cooling, also intensive massage and friction.

There have been no investigations done on the results of HYALREPAIR®-06 Bioreparant injections to women during pregnancy and lactation and children.

Administration method

The patient and the physician sign an informed voluntary consent prior to the procedure. The gel shall be injected intradermally and hypodermally. Needles to be used: 0.3x4 mm; 0.3x13 mm (30G); 0.26x4 mm; 0.26x12 mm (31G).

Recommended course:

3-6 procedures once per 2 weeks. Then one procedure per month to achieve and maintain the optimal result.

Side effects

Usage of HYALREPAIR®-06 Bioreparant can cause reactions typical for all invasive procedures: erythema, edema, slight bleeding, itch, tenderness or change of the skin colour in the place of injection. These phenomena are common after injection procedures and subside by themselves.

Interaction with other medical devices

Introduction of HYALREPAIR®-06 Bioreparant in combination with other medical devices has not been researched.

Disposal

Used gel should be disposed in compliance with local regulations.

The vial, the syringe, the needle and the remaining unused gel should be disposed of immediately after the procedure.

How supplied

5.0 ml of HYALREPAIR®-06 Bioreparant is supplied in a dark glass vial. The content is sterile. A glass vial with 2 adhesive labels for physician attached is packaged with instruction for use in a cardboard box. Apply adhesive label for physician to a patient's case history to ensure future gel identification.

Storage terms and shelf life

Store in a place protected from sunlight at the temperature between +5 and +15 °C. Do not freeze. Do not use after expiration of indicated shelf life.

DO NOT USE THE PRODUCT IF THE PACKAGE IS DAMAGED!

DO NOT STORE OPENED PACKAGE!

DO NOT RESTERILIZE!

DO NOT USE SYRINGE FOR MORE THAN ONE PATIENT!

Complaints and claims on quality of medical device are accepted by an authorized representative of the manufacturer.

Manufacturer

"Laboratory TOSKANI" LLC, Russia.

Symbols on the package

Caution!
Consult instructions for use

Do not reuse

Temperature limitation

Keep away from sunlight

Sterilized using steam or dry heat

Do not resterilize

Do not use if package is damaged

Use by date

Batch code

Manufacturer